

无菌超纯水

产品简介:

水(H₂O)是由氢、氧两种元素组成的物质, CAS 号为 7732-18-5, 分子量为 18.02; 去离子水是指除去了呈离子形式杂质后的纯水。

Leagene 无菌超纯水经去离子处理后, 经微孔滤膜过滤, 再经高压灭菌处理, 不含细菌, 几乎不含金属离子, 可用于细胞培养和 PCR 等试验, 亦可用于其他科研用途。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称	编号		
	CC0158	CC0158	Storage
无菌超纯水	100ml	500ml	RT
使用说明书	1 份		

操作步骤(仅供参考):

1、根据实验具体要求操作。

注意事项:

- Deionized Water(For Cell Culture)为无菌溶液, 请注意无菌操作。
- 试剂开封后请尽快使用, 以防影响后续实验效果。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

有效期: 12 个月有效。

相关产品:

产品编号	产品名称
CC0005	磷酸缓冲盐溶液(1×PBS,无钙镁)
CZ0063	改良台氏液(Tyrode's solution)
DH0006	苏木素伊红(HE)染色液(醇溶)
NR0003	Lezol(总 RNA 提取试剂)
PW0053	Western 抗体洗脱液(碱性)
TC0713	葡萄糖检测试剂盒(GOD-POD 比色法)

文献引用:

- Wenjuan Liu,Xiao Chen,Xiaoyong Ding,et al.Visible-light-driven cuprous oxide nanomotors with surface-heterojunction-induced propulsion.Nanoscale Horizons.January 2021.10.1039/D0NH00663G.(IF 10.989)
- Wang Yong-cai,Yang Xia,Xiao Juan,et al.Determination of the median lethal dose of zinc gluconate in mice and safety evaluation.BMC Pharmacology & Toxicology.February 2024.10.1186/s40360-024-00736-8.(IF 2.8)

注: 更多使用本产品的文献请参考产品网页