

蔗糖溶液(30%)

产品简介:

蔗糖化学式为 $C_{12}H_{22}O_{11}$, 分子量为 343.3, CAS 号为 57-50-1, 是光合作用的产物, 可作为食品添加剂, Leagene 蔗糖溶液(30%)主要由蔗糖、去离子水组成, 可用于组织脱水或其他用途, 是一种非常重要的辅助试剂。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称	编号	R00613	Storage
蔗糖溶液(30%)		100ml	4°C
使用说明书		1 份	

操作步骤(仅供参考):

- 1、根据实验具体要求操作。

注意事项:

- 1、密闭保存, 注意避免微生物污染。
- 2、为了您的安全和健康, 请穿实验服并戴一次性手套操作。
- 3、试剂开封后请尽快使用, 以防影响后续实验效果。

有效期: 4 个月有效。低温运输, 4°C 保存。

相关产品:

产品编号	产品名称
CA0075	青霉素-链霉素混合溶液(100×双抗)
DC0032	Masson 三色染色液
DA0065	台盼蓝染色液(0.4%)
DF0111	组织固定液(10% NBF)
DG0005	糖原 PAS 染色液
PS0013	RIPA 裂解液(强)
TC0699	植物总糖和还原糖检测试剂盒(DNS 比色法)

文献引用:

- 1、 Xiaoshan Liang,Xiaotao Liang,Yunyan Zhao,et al.Dysregulation of the Suprachiasmatic Nucleus Disturbs the Circadian Rhythm and Aggravates Epileptic Seizures by Inducing Hippocampal GABAergic Dysfunction in C57BL/6 Mice.JOURNAL OF PINEAL RESEARCH.July 2024.10.1111/jpi.12993.(IF 8.3)
- 2、 Wang Ruihan,Yang Yingqiu,You Guixuan,et al.Chronic lateral ankle instability using anterior tibiofibular ligament distal fascicle transfer augmentation repair: an anatomical, biomechanical, and histological study. Frontiers in Bioengineering and Biotechnology.March 2024.10.3389/fbioe.2024.1326036.(IF 4.3)

注: 更多使用本产品的文献请参考产品网页